



Drillsaw Highspeed 200™ Power System AR-200

DrillSaw Highspeed 200™ Power System

The versatile Arthrex DrillSaw Highspeed 200 power system is designed for multiple distal extremity procedures.

The unique corded power system functions at high torque and low speed, which makes it optimal for repetitive bony cuts while sparing the surrounding soft tissue. As such, this system is perfect for minimally invasive foot surgery.

The system is also compatible with all AR-300 attachments, making it a capable, corded power unit.



DrillSaw Highspeed 200™ Power System

Number One in Operation

The following can be set individually:

Program [P]

Motor speed [rpm]

Torque [Ncm]

Coolant supply volume []

Transmission ratio []



DrillSaw Highspeed 200™ Power System

Number One in Operation

Automatic
detection of motor

Ideal for minimally
invasive surgery

AR-200M can
be used with all
DrillSaw Mini
300™ attachments

Motor with cable
and motor holder
are sterilizable
(135 °C)

AR-200 Number One in the Operating Room

Powerful Motor

Top quality and power – for all surgical applications

AR-200M

Motor speed of 300 - 15 000 rpm

AR-200M-ISO

Motor speed of 300 - 40 000 rpm

Shortest and lightest 40 000 rpm motor on the market

Length: 82.5 mm

Weight: 137 g


Motor with cable and motor holder are sterilizable
(135°C)



AR-200M

Motor speed of 300 - 15 000 rpm

Coupling for AR-300 attachments




AR-200 Display

Logical Integration for Ease of Use:

One operating stage for setting the most crucial parameters

Six programs can be set individually



AR-200 Foot Pedal

Intelligent Foot Control:

Program selection



Connection and disconnection of coolant pump



Forward/reverse operation




Infinitely variable control of the motor up to the preselected speed

On/off function


IPX8- and AP-tested



AR-200 Overview - Front



AR-200 Overview - Back



AR-200 Operational Design

Operated by four buttons on the device:



Two buttons for selecting parameters




Two buttons for changing parameters



AR-200 Operational Design

Change program:




AR-200 Operational Design

Change selection:



AR-200 Operational Design

Change speed:




AR-200 Operational Design

Change selection:



AR-200 Operational Design

Torque setting:



AR-200 Operational Design

Change selection:



AR-200 Operational Design

Coolant supply volume: _____




AR-200 Operational Design

Change selection:



AR-200 Operational Design

Transmission ratio:



AR-200 Operational Design

Setup settings:

Selection menu



AR-200 Operational Design


Setup settings:

Selection menu



AR-200 Operational Design

Setup settings:



AR-200 Operational Design


Setup settings:

Selection menu




AR-200 Operational Design

Setup settings:




AR-200 Operational Design

Setup settings:



AR-200 Operational Design

Setup settings:




How to Set Up the Console



Basic Setup Procedure for AR-200

Set up the AR-200





- Connect the power cord
- Connect the foot control cable to the console



Connect the handpiece motor to the console

Basic Setup Procedure for AR-200

Irrigation Pump Tubing Setup



a
b,c
d

Open the pump arm door
Insert the irrigation pump tubing
Close the pump arm door


Basic Setup Procedure for AR-200

Follow the same sequence when removing the irrigation tubing



Connecting Attachments to the AR-200M Motor

To insert the attachment




Turn the motor chuck mechanism.




Push the adapter into the locking mechanism and release the chuck mechanism.

Connecting Attachments to the AR-200M Motor

To Remove the Attachment



Turn the motor chuck mechanism.



Pull out the adapter.


AR-200SP Spray Clip



Attach the spray clip to the motor.
The spray clip groove matches with the
coupling ring of the motor.



Insert the tubing end onto the
proximal rod end of the spray clip.



Fix the irrigation tubing in
the clips of the motor cable.

AR-200 Factory Settings

AR-200M



Program	Speed [rpm]	Torque [Ncm]	Coolant Supply Volume	Coolant	Transmission Ratio	Motor Direction
P1	15 000	100%	100%	Disabled	1:1	Clockwise
P2	8 000	100%	100%	Disabled	1:1	Clockwise
P3	3 000	100%	100%	Disabled	1:1	Clockwise
P4	15 000	100%	100%	Disabled	1:1	Clockwise
P5	15 000	100%	100%	Disabled	1:1	Clockwise
P6	15 000	100%	100%	Disabled	1:1	Clockwise

Technical Data

Console AR-200C

Power supply	110 - 130 V / 220 - 240 V, 50/60 Hz, 0.1 - 1.8 A / 0.1 - 0.9 A
Dimensions (W/H/D)	256 mm (10 in) x 109 mm (4.3 in) x 305 mm (12 in)
Weight	7.0 kg (15.4 lb)
Maximum flow rate	≥ 90 mL/min
Mains fuse	2 × 250 V – T1.6AH

AR-200M

Torque	7 Ncm
Speed	300 – 15 000 rpm



AR-200 Error Message



Button pressed during start-up



Foot control not recognized



Foot control faulty



Electronics temperature



Motor not recognized



Motor error



Motor faulty



Electronics overloaded



Service

Saw Attachments

Saw attachment – sagittal

AR-300SAG

**Saw attachment –
reciprocating**

AR-300SR



Drill Attachments

Drill attachment – style AO	AR-300DAO-2
Drill keyless chuck (0 - 3 mm)	AR-300DK30
Drill keyless chuck (2.0 - 4.5 mm)	AR-300DK45
Drill Jacobs chuck (0 - 4 mm) 1300 rpm	AR-300DJ
Drill Jacobs chuck hybrid (0 - 5 mm) 1300 rpm	AR-300DJH



Burr Attachments

Burr attachment 2.35 mm	AR-300B 15 000 rpm
Drill attachment 2.35 mm	AR-300B-2 2 750 rpm



Reamer Attachments

**Reamer adapter AO
500 rpm**

AR-300RAO

**Reamer Jacobs chuck
(0 - 4 mm) 500 rpm**

AR-300RJ



Please use the ream attachments
for hindfoot and ankle procedures!
→ Higher torque

Wire/Pin Driver Attachments

**Wire driver attachment
(0.6 - 1.6 mm)**

AR-300WD16

**Pin driver attachment
(1.0 - 2.4 mm)**

AR-300PD24



Adapters

Offset adapter

AR-200CA

Spray clip

AR-200SP



Coupling Geometries

All drill and screw attachments are offered with different coupling geometries



Round shank



Dental



Twist drill AO (mini quick)



J-notch



AR-300B

AR-300B-2



Hex (-agonal) coupling (mainly US)



AO-ASIF quick coupling (Synthes)



Color Coding

Drill attachment:










Marked green + AR-300D
1300 rpm

Ream attachment:

Marked red + AR-300R
500 rpm



Arthrex Power Family

	Prosthetics	Trauma Heavy Duty	Trauma Light Duty	Foot	Hand
<div style="background-color: #0070C0; color: white; padding: 5px; text-align: center; margin-bottom: 5px;">Possible</div> <div style="background-color: #AEC6E0; color: black; padding: 5px; text-align: center;">Optimal</div>					
 <p>DrillSaw Max 600™ power system</p>					
 <p>DrillSaw Sports 400™ power system</p>					
 <p>DrillSaw Mini 300™ power system</p>					
 <p>DrillSaw Highspeed 200™ power system</p>					



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.