


## Real-world Experience With a Decellularized Dehydrated Human Amniotic Membrane Allograft: A Prospective, Observational, Multicenter Study of a Broad Patient Population in All Wound Types. *Wounds* 2015;27(6):158-169 Janice M. Smiell, MD<sup>a</sup>; Terry Treadwell, MD<sup>b</sup>; Helen D. Hahn, RN, MBA<sup>a</sup>; Michel H. Hermans, MD<sup>c</sup>


The aim of this observational study was to gain experience in the use and performance of BIOVANCE® versus standard of care (SOC) in a real-world wound population. **A broad range of partial and full thickness wounds were studied across a mix of patient types.**

- Eligibility for inclusion included any patient that would benefit from treatment
- Unlike other chronic wound prospective, randomized, controlled trials, there were no limits on patients' age, baseline wound size or co-existing conditions
- Key comorbidities included: arterial insufficiency, autoimmune disease, diabetes, and edema/lymphedema

### BIOVANCE SUPPORTS WOUND CLOSURE ACROSS A VARIETY OF WOUND AND PATIENT TYPES



### PROJECTED WOUND CLOSURE AT 12-20 WEEKS



<sup>a</sup> Alliqua Biomedical, Inc., Langhorne, PA

<sup>b</sup> Institute for Advanced Wound Care, Montgomery, AL

<sup>c</sup> Hermans Consulting Inc, Newtown, PA

## CASE STUDY 1: VENOUS STASIS ULCER

Patient	Comorbidities	Baseline Wound Size	Closure
<ul style="list-style-type: none"> <li>61 year old female with 3 prior treatment failures</li> </ul>	<ul style="list-style-type: none"> <li>Peripheral vascular disease</li> <li>Venous insufficiency</li> <li>Immuno deficiency</li> </ul>	1.8cm x 1.2cm x 0.2cm	7 weeks with 1 application



BASELINE



WEEK 7 CLOSED

### Treatment Regimen:


- Burow's Solution was applied to the periwound area twice a day

## CASE STUDY 2: ACUTE WOUND

Patient	Comorbidities	Baseline Wound Size	Closure
<ul style="list-style-type: none"> <li>67 year old male with right BKA stump dehiscence</li> </ul>	<ul style="list-style-type: none"> <li>Severe peripheral artery disease</li> <li>Diabetes</li> </ul>	7cm x 2.5cm x 0.1cm	6 weeks with 1 application



BASELINE



WEEK 6 CLOSED

### Treatment Regimen:

- 1 application of BIOVANCE, not fenestrated
- Secondary dressing with nonadherent petroleum gauze and moist gauze

## CASE STUDY 3: DIABETIC FOOT ULCER

Patient	Comorbidities	Baseline Wound Size	Closure
<ul style="list-style-type: none"> <li>68 year old male with full thickness wound</li> </ul>	<ul style="list-style-type: none"> <li>Type 2 DM, Chronic Renal Failure, Neuropathic, Lymphedema</li> </ul>	12.9cm x 4.8cm x 0.1cm	25 weeks with 5 applications



BASELINE



WEEK 25 CLOSED

### Treatment Regimen:

- Application of BIOVANCE on Days 1, 4, 16, 38 and 136
- Secondary dressing during treatment with petroleum gauze, topical gentamicin, silver hydrofiber dressing, polyurethane foam with gauze and elastic wrap
- Patient placed on gentamicin at week 3 with positive cultures for *staphylococcus aureus*

**CITATION:** Smiell JM, Treadwell T, Hahn HD, Hermans MH. Real-world Experience With a Decellularized Dehydrated Human Amniotic Membrane Allograft. *Wounds a Compend Clin Res Pract.* 2015;27(6):158-169. <http://www.ncbi.nlm.nih.gov/pubmed/26061491>.

Presentation: The Progenerative Power of Amnion: The Science and the Clinical Experience for BIOVANCE® Human Amniotic Membrane Allograft, Mohit Bhatia, PhD. This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.

BIOVANCE is a registered trademark of Celularity, Inc. Please contact Celularity Medical Affairs to obtain additional information regarding this clinical summary at [medicalaffairs@celularity.com](mailto:medicalaffairs@celularity.com). For product information or adverse reaction reporting, telephone 1-844-963-2273. Please refer to the BIOVANCE Package Insert for complete product information. © 2021 Celularity, Inc. All rights reserved. BIO-2021-010 © 2021 Arthrex Inc. All rights reserved. DOC1-000546-en-US\_A



Marketed by Arthrex, Inc.

