

Knotless 1.8 FiberTak[®] Soft Anchor for Glenoid Labrum Repair

Surgical Technique




Knotless 1.8 FiberTak® Soft Anchor for Glenoid Labrum Repair


This tensionable knotless suture anchor combines the benefits of soft anchors with knotless soft-tissue fixation. Use the curved spear and 1.8 mm drill to precisely create a pilot hole on the glenoid rim. Insert the suture anchor through the spear, maintaining the same portal and drill trajectory. Once the suture is passed and shuttled into the locking mechanism, tension can be controlled and adjusted under direct visualization.

Advantages

- 52 lb of secure, low-profile knotless suture fixation¹
- No risk of knot impingement or loosening
- 1.8 mm drill to minimize bone removal
- Premium instrumentation with additional stability, available in curved and straight options for full access around the glenoid.
- Simple, reproducible insertion and passing techniques similar to knot-tying anchors
- Tension and retension until repair is complete



Knotless Simple Stitch




Knotless Mattress Stitch

Infinite Tissue-Tensioning Variability



Knotless 1.8 FiberTak Soft Anchor With Self-Locking Technology




Just pass it, cinch it, cut it

Knotless 1.8 FiberTak® Soft Anchor

The Knotless 1.8 FiberTak® Soft Anchor is the latest in knotless technology for instability. This design combines trusted knotless technology with the latest in suture innovation. The implant and instrument design is optimized for consistent and reliable fixation.


Implant Features



- The implant handle features a centering, stability sleeve, which increases insertion consistency by eliminating extra movement that can occur during impaction.



- The vivid, all-blue repair suture has a tapered tail.



Instrument Features





- Curved guide handles have a scallop at the proximal end for identifying the convex curve direction of the distal tip

- Flexible drill has a tri-flat design for drill hub connection, a positive stop for consistent drill depth, and a centering stability sleeve for accuracy while using the drill guide




- Laser line on the distal tip of the spear indicates the convex curve direction






Mobilize the labrum and create a bleeding bed to enhance tissue healing to bone. Pass the spear and place it on the glenoid rim. Fully advance the drill through the spear until the drill laser line or collar contacts the spear's handle.

Note: Drilling in very hard bone may require cycling the drill while maintaining consistent alignment of the drill guide. Increased size, hard bone drills are also available for use.



Insert the anchor through the spear and into bone by gentle impaction until the inserter handle is flush with the back of the spear.

Note: If insertion resistance is encountered, do not impact harder. Replace the implant and repeat the drilling/insertion process. Avoid excessive impaction as this could lead to inserter damage and/or breakage. See *warning note* on back page for additional information.




Remove the inserter handle and spear, then pull on all 3 suture tails to confirm the anchor is set in the cortical bone.


Note: A slow, steady pull is recommended to allow the anchor to properly deploy. A fast, aggressive pull could lead to improperly setting the anchor.



Retrieve the blue repair suture through the anterosuperior portal using a suture retriever. Insert a curved SutureLasso™ suture passer into the anteroinferior cannula and pass it through the capsulolabral tissue inferior to the anchor. Advance the nitinol wire loop into the joint. Retrieve the wire loop through the anterosuperior portal using the KingFisher® retriever.




5 Load the blue repair suture tail through the nitinol wire loop. Retract the wire loop through the SutureLasso™ suture passer to pull the suture to the distal end of the suture passer inside the joint. Remove the suture passer and wire loop together to shuttle the repair suture through the labral tissue.




6 Retrieve the blue repair suture and round-looped side of the white/black shuttle suture through the anterosuperior portal.

Load the repair suture through the loop of the shuttle suture. **Fold the repair suture tail at the purple mark** and crease the suture with your fingers.



7 Pull the SutureTape side of the white/black shuttle suture to transfer the repair suture back into the anchor **through the same portal where it was inserted.** **Advance the shuttle suture with repeated light tugs** until the repair suture is passed through the suture splice locking mechanism and back out of the cannula.



8 Pull the free end of the repair suture until the desired repair tension is achieved. A tissue grasper can be used to position the labrum in its desired location while applying tension on the repair. Cut the suture flush using a mini suture cutter.

Ordering Information

Implant

| Product Description | Item Number |
|-----------------------------------------------|-------------|
| Knotless 1.8 FiberTak® Soft Anchor, #2 Suture | AR-3636 |



Disposable Instruments

| Product Description | Item Number |
|-------------------------------------------------------------------------------------------------------|--------------|
| Percutaneous Instrument Kit for Knotless FiberTak Soft Anchor, w/ 1.8 mm rigid drill | AR-3610PK-3 |
| Knotless FiberTak Disposable Kit, w/ tapered curved spear, 1.8 mm flexible drill, and blunt obturator | AR-3610DC-3 |
| Knotless FiberTak Disposable Kit, w/ curved spear, 1.8 mm flexible drill, and blunt obturator | AR-3638DC |
| Knotless FiberTak Disposable Kit, w/ straight spear, 1.8 mm rigid drill, and blunt obturator | AR-3638DS |
| 1.8 mm Flexible drill w/ hub, trocar obturator, sterile | AR-3610ND-2 |
| 1.9 mm Flexible drill w/ hub, trocar obturator, sterile, hard bone | AR-3610ND-4 |
| 1.8 mm rigid drill, sterile | AR-3600D-2 |
| 1.8 mm flexible ShaverDrill™ device | AR-3610NSD-2 |

Reusable Curved Spear Options

| Product Description | Item Number |
|------------------------------------------------------------------|--------------|
| Curved Spear for FiberTak soft anchor, w/ trocar obturator | AR-3610CTC |
| Tight Curved Spear for FiberTak soft anchor, w/ trocar obturator | AR-3610CTC-2 |

Reusable Straight Spear Options

| Product Description | Item Number |
|--------------------------------|-------------|
| Fishmouth spear | AR-3610F |
| Spear w/ circumferential teeth | AR-3610CT |
| Slotted spear | AR-3610ST |



Fishmouth Spear
AR-3610F



Spear w/
Circumferential Teeth
AR-3610CT



Slotted Spear
AR-3610ST

Ordering Information (Cont)

Knotless 2.6 FiberTak® Soft Anchor

With the same tensionable knotless technology as the smaller Knotless 1.8 FiberTak soft anchor, the Knotless 2.6 FiberTak soft anchor includes a #5 repair suture. This anchor option combines the benefits of a soft anchor with a broader soft-tissue repair option. Using a drill guide and 2.6 mm drill, create a pilot hole and insert the anchor through the drill guide. Once the repair suture is passed through tissue, shuttle it into the knotless suture mechanism. Suture repair tension can be controlled and adjusted under direct visualization.



Implant

| Product Description | Item Number |
|-------------------------------------------------|-------------|
| Knotless 2.6 FiberTak soft anchor, w/ #5 suture | AR-3641 |

Disposable Instruments


| Product Description | Item Number |
|--------------------------------------------------------|-------------|
| Disposables Kit, w/ spear, obturator, and 2.6 mm drill | AR-3650DS |
| 2.6 mm drill | AR-3657 |
| 2.6 mm ShaverDrill™ device | AR-3657SD |

Reusable Instruments

| Product Description | Item Number |
|-------------------------------------------------------|-------------|
| Punch for Knotless 2.6 FiberTak soft anchor | AR-3656 |
| Blunt-Tip Obturator for Spear | AR-3658B |
| Trocar-Tip Obturator for Spear | AR-3658T |
| Angled Spear w/ circumferential teeth | AR-3655 |
| Circumferential teeth Spear, 2.6 FiberTak soft anchor | AR-1941CT |
| Fishmouth Spear, 2.6 FiberTak soft anchor | AR-1941DGF |



Punch for Knotless 2.6 FiberTak Soft Anchor
AR-3656




Blunt-Tip Obturator
AR-3658B




Trocar-Tip Obturator
AR-3658T



Angled Spear w/
Circumferential Teeth
AR-3655



Spear w/
Circumferential Teeth
AR-1941CT



Fishmouth Spear
AR-1941DGF

Reference

1. Arthrex, Inc. Data on file (APT 3531). Naples, FL; 2017.

WARNING!

TO HELP AVOID INSERTER BREAKAGE AND POTENTIAL PATIENT INJURY:

- Avoid excessive impactation as this could lead to inserter damage and/or breakage.
- If insertion resistance is encountered, do not impact harder. Replace the implant and repeat the drilling/insertion process.
- Visually inspect the inserter for potential breakage after each implantation. See image below for reference **(a)**.



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.



Arthrex Manufacturer,
Authorized Representative
and Importer information:
eIFU - Arthrex



View U.S. patent information

arthrex.com