


# Spare and Replacement Parts for Synergy Scopes

Quick Guide

AR-3350\* | AR-3351\* | AR-3352\* | AR-3355\*



O-ring kit for  
C-mount scopes (5 per kit)  
**8300020859**



Wolf adapter  
**8300042617**



Storz adapter  
**89542801**



Polishing paste  
**8101001000**

#### Attention:

Polishing paste can be used on the proximal and distal glass surfaces of the scopes to aid in removing debris or small scratches.

- If the scope is removed from body, the light must be placed on standby.
- The scope and light cord should never be placed on anything flammable unless the light is on standby.
- If you are experiencing a hot light post junction, please inspect the light cords and light post for debris or damage.

\*For further information on our item numbers, please refer to the product catalog.

This is not medical advice and Arthrex recommends that surgeons be trained in the use of any particular product before using it in surgery. A surgeon must always rely on their own professional judgment when deciding whether to use a particular product when treating a particular patient. A surgeon must always refer to the package insert, product label, and/or directions for use before using any Arthrex product. Products may not be available in all markets because product availability is subject to the regulatory and medical practices in individual markets. Please contact your Arthrex representative if you have questions about the availability of products in your area.

[arthrex.com](http://arthrex.com)

© 2022 Arthrex, Inc. All rights reserved. LI1-00045-EN\_B

