


Nano Arthroscopic Instrumentation



Arthrex® 



Nano Instrumentation

Size Matters in Arthroscopic Resection

Introducing Nano single-use, sterile-packed arthroscopy instrumentation.

Harnessing 20 years of excellence in designing arthroscopic hand instrumentation, Arthrex has engineered the next generation of instruments for diagnostic, resection, and extraction procedures.

The NanoProbe is an articulating, retractable hook probe used for diagnostic procedures. The NanoBiter, NanoScissor, and NanoGrasper are small arthroscopic soft-tissue resection instruments. Just 2 mm in diameter, they are designed for atraumatic insertion through tight joint spaces. These instruments are sharp and engineered to efficiently resect meniscal tissue and other soft-tissue structures.

Key Features and Benefits

■ 2 mm, nano-sized instrumentation

- Designed for use in traditional arthroscopic operative theaters or procedure rooms, or to complement in-office arthroscopy.
- Can be used with a 4 mm arthroscope, small-joint arthroscope, or 1.9 mm NanoScope™ system.
- Available in standard (130 mm) or small-joint (70 mm) working lengths.

■ More effective, atraumatic resection


- Low-profile design enables greater access to hard-to-reach pathology without sacrificing resection efficiency.
- Convenient, single-use design allows for an unmatched cutting-edge sharpness that ensures optimal resection on every procedure.
- Provided sterile for a precise, clean instrument for every patient.

Ordering Information


Knee, Shoulder, and Hip Procedures	Item Number
NanoProbe, disposable (single-pack)	AR-10100N
NanoGrasper, straight, 130 mm, disposable (single-pack)	AR-10913D-1
NanoScissor, straight, 130 mm, disposable (single-pack)	AR-10915D-1
NanoBiter, straight, 130 mm, disposable (single-pack)	AR-10911D-1
NanoBiter, 15° up, 130 mm, disposable (single-pack)	AR-10922D-1

Small-Joint Procedures	Item Number
NanoProbe, disposable (single-pack)	AR-10100N
Small-Joint Probe, reusable	AR-30000
NanoBiter, straight, 70 mm, disposable (single-pack)	AR-10901D-1
NanoBiter, 15° up, 70 mm, disposable (single-pack)	AR-10902D-1
NanoGrasper, straight, 70 mm, disposable (single-pack)	AR-10903D-1
NanoScissor, straight, 70 mm, disposable (single-pack)	AR-10905D-1


NanoBiter



3.4 mm Straight Punch



Bite Size Comparison



For more information about Nano Arthroscopic Instrumentation or the NanoScope Arthroscopy System, please contact your local Arthrex representative or visit Arthrex.com.

www.arthrex.com

© 2019 Arthrex, Inc. All rights reserved. FL1-000036-en-US_C

