


# Univers Apex OptiFit™ Total Shoulder System

Optimized Stemmed Anatomic Arthroplasty



- In situ adjustability with variable inclination from 125°-140° and version  $\pm 10^\circ$
- Rectangular stem geometry provides rotational stability
- Titanium plasma spray enables ingrowth fixation
- 85 mm effective length optimizes stability while sacrificing minimal bone
- Optimized subscapularis and revision techniques



Learn more about the  
OptiFit total shoulder system



© 2024-12 Arthrex, Inc. All rights reserved.  
FL1-000930-en-US\_A