



# Fluid Containment

Keep Your OR Dry and Slip Free

## EcoDri-Safe™\* Absorbent Rolls



Containing fluids and streamlining cleanup procedures has never been easier than with EcoDri-Safe absorbent rolls. The 28 in × 100 ft rolls are perforated every 20 inches to create 60 easy-to-tear-off 20 in × 28 in floor mats on every roll.

Product Description	Item Number
EcoDri-Safe Rolls, non-poly-backed, 100 ft × 28 in, 20 in × 28 in, 60 sheets per roll	USA-ECO100
EcoDri-Safe Rolls, poly-backed, 100 ft × 28 in, 20 in × 28 in, 60 sheets per roll	USA-ECO100PLY
EcoDri-Safe Dispensing Cart	USA-ECODISP

## The Suctioner®\* Mat



This antimicrobial surgical mat incorporates powerful, quiet suction capability with exceptional anti-fatigue properties of The Surgical Mat™ floor covering and the absorbency of the Dri-Safe® pad.

Product Description	Item Number
The Suctioner Surgical Mat, 22 in × 35 in (comfort mat), 42 in × 56 in (total), blue, qty. 8	USA-X10

## The EcoSuctioner®\* Mat



- All-in-one ergonomic suction mat that allows two people to stand comfortably during a procedure
- Designed with directional grooves for improved quiet, continual suction; attaches to any suction source
- Large, rugged absorbent mat included on underside creates more nonslip surface in OR

Product Description	Item Number
The EcoSuctioner Surgical Mat, 28 in × 40 in (comfort mat), 42 in × 56 in (total), green, qty. 8	USA-XG8

\*EcoDri-Safe, The Surgical Mat, Dri-Safe, The SuctionBlu, The Suctioner, The SuctionBlu, and The EcoSuctioner are trademarks of Flagship Surgical, LLC.

## The SuctionBlu™\* Mat



SuctionBlu mat is a smaller, more economical version of The EcoSuctioner mat but larger than The Suctioner mat.

Product Description	Item Number
The SuctionBlu Surgical Mat, 42 in × 56 in, blue	USA-XBL



## Fluid Diamond®\* Floor Aspirator



- Fast, effective fluid collection
- Powerful suction action keeps floors dry and reduces the potential for slipping or falling in the OR
- Attaches to any surgical suction source
- Includes suction tubing for easy attachment

Product Description	Item Number
Fluid Diamond Floor Aspirator, qty. 30	USA-FDX

\*Fluid Diamond is a trademark of Flagship Surgical, LLC.


## Dri-Safe® Absorbent Pads



- Superabsorbent, economical, disposable floor pads that eliminate expensive laundry costs
- Dri-Safe pads lock-in fluids to support a safer OR and are available in poly- and non-poly-backed options

Product Description	Item Number
Dri-Safe Absorbent Pad, 28 in × 40 in, blue, qty. 20	USA-B28X40
Dri-Safe Absorbent Pad, 28 in × 40 in, blue, poly-backing, qty. 20	USA-B28X40PLY
Dri-Safe Absorbent Pad, 28 in × 56 in, blue, qty. 15	USA-B28X56
Dri-Safe Absorbent Pad, 28 in × 56 in, blue, poly-backing, qty. 15	USA-B28X56PLY

## The Surgical Mat™ and The Mini Mat® Surgical Mats



- The Surgical Mat and The Mini Mat surgical mats have resilient memory foam for enhanced comfort
- Nonslip design to promote a safer OR
- Custom tapered fluid channels to help keep the surgical team's feet dry and reduce risk of cross contamination

Product Description	Item Number
The Surgical Mat Floor Covering, 20 in × 39 in, blue, qty. 10 (a)	USA-010
The Mini Mat Surgical Mat, 13 in × 17 in, blue, qty. 10 (b)	USA-A10