

Arthrex Screwdriver Reference Chart

for TITANIUM INTERFERENCE SCREWS

	I399	I586R	I896	I995GT	I995SHN	I998	I998GT	I998SH	I998T
5 mm Interference Screws (noncannulated)					X				
6 mm Interference Screws				X	X				
7, 8, 9 & 10 mm Interference Screws[^] and Soft Screws[^]	X		X			X	X	X	X
Titanium Tibial and Femoral RetroScrews (cannulated)		X							
Low Profile Cancellous Screws (noncannulated) and Bi-Cortical Posts (noncannulated)					X				

 **AR-I399** Pinlock I Cannulated Screwdriver, 3.5 mm Hex[^]

 **AR-I586R** RetroScrew Driver, thin

 **AR-I896** Pinlock II Cannulated Screwdriver, 3.5 mm Hex[^]

AR-I995GT* Guide Pin Tip Transtibial Screwdriver Shaft, 2.5 mm Hex


AR-I995SHN* Noncannulated, Short Screwdriver Shaft, 2.5 mm Hex

AR-I998* Cannulated Screwdriver Shaft, 3.5 mm Hex[^]

AR-I998GT* Guide Pin Tip Transtibial Screwdriver Shaft, 3.5 mm Hex

AR-I998SH* Cannulated, Short Screwdriver Shaft, 3.5 mm Hex

AR-I998T* Cannulated Transtibial Screwdriver Shaft, 3.5 mm Hex

* Requires Ratcheting Screwdriver Handle (AR-I999), Non-Ratcheting Screwdriver Handle (AR-I999NR) or Tear Drop Handle (AR-2001) 

[^] Uses 2 mm Nitinol Guide Pin (AR-I254)



Arthrex Screwdriver Reference Chart

for BIOABSORBABLE INTERFERENCE SCREWS & POSTS

	I357DB	I367DB	I386	I486	I586R	I997	I997D	I997DSH	I997GT	I997SHG	I997T	I996CD	I996CD-I
Bio-Interference Screws^{^^} (cannulated)			X			X			X	X	X		
Full Thread Bio-Interference Screws^{^^}			X			X					X		
Delta Bio-Interference Screws^{^^} (cannulated), Proximal Bio-Cortical Screws^{^^} (cannulated) and Distal Bio-Cortical Screws (noncannulated)			X	X		X	X	X			X		
Tibial and Femoral RetroScrews (cannulated)					X								
Bi-Cortical Bio-Posts (noncannulated)		X											
Bio-Post & Washer (noncannulated)	X												
BioComposite Interference Screws (full thread, round delta tapered & delta tapered)												X	X

AR-I357DB* Bio-Post Driver

AR-I367DB* Bi-Cortical Bio-Post Driver

 **AR-I386** Cannulated Screwdriver for Bio-Interference Screw^{^^}

 **AR-I486** Cannulated Screwdriver for Delta Bio-Interference Screw^{^^}

 **AR-I586R** RetroScrew Driver, thin

 **AR-I996CD** BioComposite Interference Screw Driver

 **AR-I996CD-I** BioComposite Interference Screw Driver, quick connect

AR-I997* Cannulated Screwdriver Shaft for Bio-Interference Screw


AR-I997D* Cannulated Screwdriver Shaft for Delta Bio-Interference Screw

AR-I997DSH Cannulated Short Screwdriver Shaft for Delta Bio-Interference Screw^{^^}


AR-I997GT* Guide Pin Tip Transtibial Screwdriver Shaft for Bio-Interference Screw

AR-I997SHG* Guide Pin Tip Short Screwdriver Shaft for Bio-Interference Screw

AR-I997T* Cannulated Transtibial Screwdriver Shaft for Bio-Interference Screw^{^^}

* Requires Ratcheting Screwdriver Handle (AR-1999), Non-Ratcheting Screwdriver Handle (AR-1999NR) or Tear Drop Handle (AR-2001) 

^{^^} Uses 1.0 mm Nitinol Guide Pin (AR-1249)

 Color bar references the color of the screwdriver handle