


Calcaneal Fractures

Innovative Treatment Solutions


Perimeter Plates

- Low-profile anatomic locking plates, 1.35 mm thick
- Commercially pure titanium for easy contouring
- 3.5 mm nonlocking and variable-angle locking screws
- Right and left plates in XS, S, M, and L sizes





Percutaneous Plates

- Five low-profile locking plate designs
- Attachable percutaneous guides
- Tapered ends to aid in insertion
- Right and left plates with standard and long sizes




5.0 mm/7.0 mm Compression FT Screws

- Cannulated, fully threaded headless compression screws
- Variable-stepped thread pitch provides continuous compression
- All instruments and implants in a single tray



Ordering Information


Product Description	Item Number
Calcaneal Fracture Implants and Guides	AR-8950S-05
Calcaneal Fracture Perimeter Plates	
X-Small, left	AR-8946-01
Small, left	AR-8946-02
Medium, left	AR-8946-03
Large, left	AR-8946-04
X-Small, right	AR-8946-05
Small, right	AR-8946-06
Medium, right	AR-8946-07
Large, right	AR-8946-08
Calcaneal Fracture Percutaneous Plates	
Anterior/Posterior, standard, left	AR-8946-01
Anterior, standard, left	AR-8946-02
Anterior/Posterior, long, left	AR-8946-03
Anterior, long, left	AR-8946-04
Anterior/Posterior, standard, right	AR-8946-05
Anterior, standard, right	AR-8946-06
Anterior/Posterior, long, right	AR-8946-07
Anterior, long, right	AR-8946-08
Low Profile T-Plate, 4-hole	AR-8954MT
Calcaneal Fracture Instruments	AR-8950S-04
5.0 mm/7.0 mm Compression FT Screw System	AR-8750S



Perimeter Plate



Percutaneous Anterior Process Plate




Percutaneous Anterior Process/Posterior Tuberosity Plate




Low Profile T-Plate

Calcaneal Fracture Instruments
AR-8950S-04



5.0 Large
Compression FT Screw



7.0 X-Large
Compression FT Screw



Quickset™ Calcium Phosphate Cement

*Quickset is a trademark of Graftys, S.A.