




SixFix Circular Frame System




Location	Description
1	SixFix Ring Case
2	SixFix Supplemental Fixation Case
3	SixFix Instrument Tray
4	SixFix Fixation Element Case
5	SixFix Foot Rockers, Short
6	4,5,6mm Half Pins

76-11607: Ring Case



Location	Part Number	Description	QTY
1	71-6120L	Foot Ring, 180mm, Long	1
1	71-6120S	Foot Ring, 180mm, Short	1
2	71-6100L	Foot Ring, 150mm, Long	1
2	71-6100S	Foot Ring, 150mm, Short	1
3	71-90150	Foot Arch, 150mm	1
3	71-90180	Foot Arch, 180mm	1
4	71-10135	Full Ring, 135mm	2
4	71-58135	2/3 Ring, 135mm	1
5	71-10155	Full Ring, 155mm	2
5	71-58155	2/3 Ring, 155mm	1
6	71-10175	Full Ring, 175mm	2
6	71-58175	2/3 Ring, 175mm	1
7	71-10195	Full Ring, 195mm	2
7	71-58195	2/3 Ring, 195mm	1


76-11603: Supplemental Fixation Element Case



Location	Part Number	Description	QTY
1	76-10104	Wire Fixation Bolt, XL	30
2	76-10149	Half Pin Fixation Bolt	10
3	76-10180	6mm Threaded Nut, 10mm Hex	60
4	76-10182	Flange Nut, 10mm Hex	60
5	76-10185	Hex Nut M6, Tall	10
6	76-10161	Hex Bolt 6 mm X 12 mm	25
7	76-10162	Hex Bolt 6 mm X 16 mm	25
8	76-10163	Hex Bolt 6 mm X 20 mm	25
9	76-10171	Flange Hex Bolt 6 mm X 12 mm	20
10	76-10172	Flange Hex Bolt 6mm X 16mm	20
11	76-10173	Flange Hex Bolt 6 mm X 20 mm	20
12	76-10203	Rancho Cube, 4-Hole	4
13	76-10204	Rancho Cube, 5-Hole	4
14	76-10125	2-Hole Threaded Posts	5
15	76-10140	4-Hole Threaded Posts	4
16	76-10145	5-Hole Threaded Posts	4
17	76-10146	6-Hole Threaded Posts	4
18	76-11085	Threaded Rods, 85mm	8
19	76-11115	Threaded Rods, 115mm	8
20	76-11200	Threaded Rods, 200mm	8
21	76-11300	Threaded Rods, 300mm	8
22	76-11450	Transfixing Bone Screw	2
23	76-10032	10mm Wrench, Low Profile	1
23	76-10018	L-Wrench, 5mm Hex w/ Ball	1
23	76-10030	L-Nut Driver, 10mm	1
23	76-10032	10mm Angle Wrench, Low Profile	1




76-11608: Fixation Element Case, Tray 1



Location	Part Number	Description	QTY
1	76-10100	Variable Angle Clamp, Short	3
1	76-10101	Variable Angle Clamp, Tall	2
2	76-10030	L-Nut Driver, 10mm	1
3	76-10020	Combo Wrench, 10 mm	2
4	76-10023	Needle Nose Pliers	1
5	76-10066	Screw Guide/Tissue Protector	1
6	76-10044	Pin Adaptor QC AO 4mm (Blue) AMDT #	1
6	76-10054	Pin Adaptor QC AO 5mm (Yellow) AMDT #	1
6	76-10064	Pin Adaptor QC AO 6mm (Black) AMDT #	1
7	76-10035	Universal Chuck w Handle	1
8	71-61150	Butt Frame Connector, 150mm	2
9	71-61180	Butt Frame Connector, 180mm	2
10	76-10191	Extended Connector, Short	8
11	76-10192	Extended Connector, Long	8
12	76-10165	Rancho Bolt 6mm x 8mm	10
13	76-10209	Centering Sleeve, 4mm	10
14	76-10210	Centering Sleeve, 5mm	10
15	76-10112	Washer	20
16	76-10203	Rancho Cube, 4-Hole	4
17	76-10204	Rancho Cube, 5-Hole	4
18	76-10207	Rancho w/Post, 2-Hole	4
19	76-10208	Rancho w/Post, 3-Hole	4




76-11608: Fixation Element Case, Tray 2




Location	Part Number	Description	QTY
1	76-10040	Fixation Pin Drill Sleeve	2
2	76-10045	Trocar	2
3	76-10043	Drill Bit, 3.0mm X 250mm	2
4	76-11506	Static Strut, Short (30mm)	3
5	76-11507	Static Strut, Medium(40mm)	3
6	76-11508	Static Strut, Long (60mm)	5
7	76-11509	Static Strut, Extra Long (80mm)	5
8	76-10000	Wire Tensioner	1
9	641018400	Smooth Wire 1.8x400mm	15
9	643518400	Olive Wire 1.8x400mm	15

76-11606: Instrument Case, Tray 1



Location	Part Number	Description	QTY
1	76-10000	Wire Tensioner w/ Nose	1
2	76-10010	Modular TL Handle, AO Connect	2
3	76-10011	Hex Driver, 3/32", AO QC	2
4	71-00002	Strut Screw Add-on, 100mm	6
5	25-0040	Shoulder Screw	3
6	76-10012	Hex Driver, 5mm, Fiducial	2
7	76-10020	Combo Wrench, 10 mm	2
8	76-10019	COMBO Wrench, 13mm	1

76-11606: Instrument Case, Tray 2



Location	Part Number	Description	QTY
1	73-00001	SixFix Strut Assembly, 50mm - Base Unit	1
2	71-00001	SixFix Strut Assembly, 100mm - Base Unit	1
3	73-00002	Strut Screw Add-on, 50mm	6



XXXX_STAR 90 F4
 External Fixator Cortical
 Bone
 Screws_240104_PK_026_
 CL.psd

Package Inventory

71-6130S_AM
 DT
 Footlockers
 and
 Clips_240104_
 TT_020_CL.psd



Location	Part Number	Description	QTY
1	71-00004	Strut Clip, Plastic	6
2	71-6130S	Foot Rocker, Short (pair)	1
3	F4-134120	Half Pin 4x34x120mm	4
3	F4-134150	Half Pin 4x40x150mm	4
3	F4-145120	Half Pin 5x34x120mm	4
3	F4-145150	Half Pin 5x40x150mm	4
3	F4-145180	Half Pin 5x50x180mm	4
3	F4-156150	Half Pin 6x40x150mm	4
3	F4-156180	Half Pin 6x50x180mm	4
3	F4-156200	Half Pin 6x60x200mm	4
3	ST134120HA	Half Pin 4x34x120mm HA Coated	4
3	ST134150HA	Half Pin 4x40x150mm HA Coated	4
3	ST145120HA	Half Pin 5x34x120mm HA Coated	4
3	ST145150HA	Half Pin 5x40x150mm HA Coated	4
3	ST145180HA	Half Pin 5x50x180mm HA Coated	4
3	ST156150HA	Half Pin 6x40x150mm HA Coated	4
3	ST156180HA	Half Pin 6x50x180mm HA Coated	4