

# The SuctionBlu™\* Mat

Combining Anti-fatigue Properties  
With Suction Capability



The SuctionBlu disposable surgical mat combines the absorbency of the poly-backed Dri-Safe®\* pad and the nonslip, anti-fatigue memory foam construction of The EcoSuctioner®\* mat with powerful, quiet, continuous suction capability. These features add up to ensure a drier, safer, and cleaner operating room.

The graduated slots on The SuctionBlu mat use gravity to direct fluids toward the suction tube. The SuctionBlu mat easily attaches to any operating room vacuum

source in seconds and efficiently manages unwanted fluids, keeping them off the OR floor.

The SuctionBlu mat is a smaller, more economical version of The EcoSuctioner mat but it is larger than The Suctioner®\* mat. The chart below compares the two suction mats by size and weight.

## Disposable Anti-fatigue Mats

Item Number	Product Description	Size of Surgical Mat	Total Size With Dri-Safe Pad	Weight
USA-XBL	The SuctionBlu Mat	25 in × 37.5 in	42 in × 56 in	1 lb 6 oz
USA-XG8	The EcoSuctioner Mat	28 in × 40 in	42 in × 56 in	1 lb 10 oz
USA-X10	The Suctioner Mat	22 in × 35 in	42 in × 56 in	4 lb 12 oz

\*The SuctionBlu, Dri-Safe, The EcoSuctioner, and The Suctioner are trademarks of Flagship Surgical, LLC.